

# Improve mRNA translation and stability with customizable poly(A) tailing

## A-Plus™ Poly(A) Polymerase Tailing Kit

### Product Intro

The A-Plus™ Poly(A) Polymerase Tailing Kit is for post-transcriptional poly(A) tail addition to RNAs, enabling *in vitro* tailing to lengths equal to or greater than those found in nature. The standard protocol produces a poly(A) tail length of ~150 bases on 60 µg of RNA. The protocol can be easily adjusted to obtain tail lengths greater than 300 bases which cannot be accomplished with co-transcriptional tailing methods.


### Benefits

- **Customizable:** Generate poly(A) tails longer than those in an encoded tail, even >300 bases.
- **Increased stability:** Longer poly(A) tails increase RNA stability and translation efficiency.
- **Efficient workflow:** Directly input capped RNA from [ScriptCap™](#) reactions with no clean up.

### Product Description

The [A-Plus™ Poly\(A\) Polymerase Tailing Kit](#) uses ATP as a substrate for template-independent, post-transcriptional addition of a poly(A) tail to the 3'-end of RNA. The standard protocol produces a poly(A) tail length of ~150-200 bases on 60 µg of RNA, but tail lengths of >300 bases are possible using a modified protocol. Polyadenylation increases the stability of mRNA in eukaryotic cells and enhances its ability to be translated after transfection or microinjection. Additionally, a poly(A) tail is useful to provide a priming site for first-strand cDNA synthesis in certain applications and can be used to end-label or quantify mRNA.

### Performance Data



**Figure 1.** Different poly(A) tail lengths generated using the A-Plus™ Poly(A) Polymerase Tailing Kit. The standard reaction results in a poly(A) tail length of approximately 150 bases, but longer tail lengths can be achieved by adjusting enzyme concentration, incubation time, RNA input, and reaction volume.

### Ordering information

| Catalog Number | Description                                           |
|----------------|-------------------------------------------------------|
| C-PAP5104H     | A-Plus™ Poly(A) Polymerase Tailing Kit (50 reactions) |

Visit [cellscript.com](https://www.cellscript.com)  
for more information



For research use only